

Restriction mapping

Aim:- To construct restriction maps using NEB cutter version 2.0

Query:- P

Introduction:-

This tool will take a DNA sequence & find the large non-overlapping open reading frames using the E. coli genetic code & sites for all Type II & commercially available type III restriction enzymes that cut sequences just once. By default only enzymes available from NEB are used. But other sets may be chosen. Just enter your sequence & submit. Further option will appear with output. The maximum size of the input is 1 Mbp & the maximum sequence length is 300 kbases.

Query:- PBR322

URL:- www.nebi.nlm.nih.gov/nc2.neb.com/ncb_cutter2/

Procedure:-

- log on to NCBI (www.ncbi.nlm.nih.gov/)
- Retrieve a nucleotide sequence of your own interest in FASTA format.

- Log on to (nc2.web.com) NEB cutter 2/)
- Paste the retrieved nucleotide seq in the space given in NEB cutter welcome page.
- Click on either linear (Circular option.
- Click on either NEB enzyme / commercial enzyme option.
- Click on Submit option.
- Phase I result will be displayed.

Procedure for phase II Results:-

- Click onto the option called 2 cutter/3 Cutter on phase I result page.
- On clicking into 2 cutters/3 Cutter the tool shows about the no. of times a restriction enzyme cuts.
- Select two points on sequence obtained on phase I results & click on flanking enzyme option.
- This phase II result are obtained.

Procedure for phase III Results:-

- Click on option called custom digest
- Then select the enzyme of your interest & click onto option called custom digest.
- In order to view gel electrophoresis which on a option 'view gel' that gives the gel electrophoresis diagram.
- In the gel electrophoresis diagrammatic page. User has an option to given their interest

Name of Experiment.....
Experiment No. 6

Date 17/1/24
Experiment Result.....

Page No. 21

of marker agarose 1%
Thus phase in results is obtained.

Result & Discussion:-

A restriction mapping is done with the help of NEB cutter V2.0 tool.

This tool not only used to obtain the restriction map but also action / the user to obtain the following data.

- i) This tool allow the user to obtain the result in both linear & Circular pattern.
- ii) 2 Cutter :- A particular restriction enzyme which has two restriction sites on DNA.
- iii) 3 Cutter: A particular restriction enzyme cleaves the DNA thrice.
- iv) Flanking enzymes.
- v) Custom digest [This option allow the user to select main interested restriction enzymes list. in order to clean the given DNA query].
- vi) View gel [This option allow the user to view electrophoresis pattern of the restriction maps which have been obtained custom digest].